

UNICSKIN

YOUR EFFECTIVE SOLUTION

UNICTECH EYE MASK

DEVICE WITH LED LIGHT & EMS TECHNOLOGIES
TRIPLE ACTION: ANTI-WRINKLE + ANTI-PUFFINESS +
ANTI-DARK CIRCLES



“UNICSKIN beats the aesthetic center with this revolutionary device for the eye contour.”

“100% ophthalmologically safe”

Dr. Sierra, n°1 in Oculoplastic Surgery in New York.

USER MANUAL



LET'S BEGIN

Congratulations! You've taken the first step towards youthful, radiant eyes by choosing UNICSKIN, global experts in techno-beauty devices that allow you to care for your skin effectively in your own home.

Before you start using this device, please read all instructions and precautions carefully and only use it for the purpose described in this manual.

1 DEVICE, 2 TECHNOLOGIES, 3 ACTIONS

UNICTECH EYE MASK is designed as a safe eye contour treatment for home use. **It combines advanced electrical muscle stimulation (EMS) technology with LED light to stimulate and activate the skin cells of the eye contour.** It improves microcirculation and promotes collagen synthesis to reduce wrinkles, puffiness and dark circles while improving skin elasticity and firmness.

TAKE A CLOSER LOOK AT UNICTECH EYE MASK

LED LIGHT TECHNOLOGY (RED LIGHT)

Long wave light (630 nanometres) that penetrates deep into the skin to slow down ageing and correct expression lines by activating collagen and elastin formation.

MICROCURRENT TECHNOLOGY (EMS)

Microcurrents mimic the body's electrical currents at cell level to tone, firm and reduce wrinkles around the eyes.

LIGHT ERGONOMIC DESIGN

Fits the shape of your face and eye contour perfectly, so you can use it while doing other things.

ON/OFF BUTTON

Press and hold the power button for one second to turn on the device.

EMS/LED INTENSITY

You can customise the intensity of the treatment. There are three different intensity levels: low, medium and high.

- **LOW INTENSITY:**
MILD WRINKLES/PUFFINESS/DARK CIRCLES
- **MEDIUM INTENSITY:**
MODERATE WRINKLES/PUFFINESS/DARK CIRCLES
- **HIGH INTENSITY:**
SEVERE WRINKLES/PUFFINESS/DARK CIRCLES

USB CHARGING PORT

One charge is enough for ten uses.



COMPLETE YOUR ROUTINE

We recommend our UNICEYES TRIPLE ACTION eye contour cream and UNICCOLLAGEN FLASH EYE MASK patches.



TIMER

Turns off after ten minutes to indicate the end of the treatment.



USB CHARGING CABLE

Charge the device anywhere, any time.



ACCESSORIES

Conveniently store the device and take it anywhere.

HOW TO USE UNICTECH EYE MASK



STEP 1

We recommend **applying your daily beauty routine** to clean skin. First, apply our **UNICEYES TRIPLE ACTION** serum to the eye contour area. The eye contour area needs to be slightly moisturised to enhance the effect of the EMS technology.



STEP 2

Remove the device from the pouch and press the power button for one second to turn it on. You'll see the red LEDs on both sides light up.



STEP 3

Adjust the treatment intensity (three different intensities) by pressing the EMS + LED button.

- **LOW INTENSITY:**
MILD WRINKLES/PUFFINESS/DARK CIRCLES
- **MEDIUM INTENSITY:**
MODERATE WRINKLES/PUFFINESS/DARK CIRCLES
- **HIGH INTENSITY:**
SEVERE WRINKLES/PUFFINESS/DARK CIRCLES



STEP 4

Place the glasses over the eye contour area for ten minutes. The glasses will turn off automatically when the time is up.

The LED light and EMS technology work upwards. The shape and design of the glasses allow the current to activate and move up both sides of the face before reaching the crow's feet, releasing the effects of the light and microcurrent at lower levels over this delicate area, particularly the eyelids.

For optimal results, we recommend applying the UNICTECH EYE MASK glasses three times a week for four weeks. You can use them anywhere, any time - whenever suits you. **The effects can be seen from the very first session.**

One month is enough to see optimal results and improvements in the eye contour.

After twelve sessions, we recommend applying UNICTECH EYE MASK glasses once or twice a week to maintain the treatment results.



UNICSKIN recommends ten-minute sessions to achieve the perfect synergy and effectiveness of LED and EMS technology in the eye area. Our technology has been tested and proven by experts and fully meets the world's highest health standards.

CHARGING

Your UNICTECH EYE MASK device comes with a USB charging cable. A one-hour charge will be enough to fully recharge the device and provide up to ten ten-minute treatments.

While recharging, the indicator lights will flash red. When the device is fully charged, the three lights will turn green.

HOW TO LOOK AFTER THE DEVICE

After use: Clean the device with the cleaning cloth provided or with a slightly damp cloth (the device must not come into contact with water. See PRECAUTIONS section).

Store the device in the velvet pouch provided and then in the hard case to keep it protected and take it anywhere.

TECHNICAL SPECIFICATIONS

Size	139.6 x 53.5 x 37mm
Weight	42g
Volts	3.7V
Frequency	0.5W
Charging volts	5V
Charging current	120mA
Max. Operating current	120mA
EMS Frequency	Continuous mode 1KHZ



PRECAUTIONS

- It is normal to feel slight electric shocks while using the device.
- Not water resistant.
- Do not wash with water.
- Store in a cool ventilated place.
- Do not dismantle the device.
- Not recommended in pregnancy.
- Keep out of reach of children.



Caution: You may feel slight painless shocks in the eye contour area during treatment. If you find this uncomfortable, use a lower intensity

WARRANTY

No refunds will be given on a used product purchased on our website (www.unicskin.com) unless due to a manufacturing fault. If the product was not purchased on our website, please refer to the retailer's warranty policy, which under no circumstances replaces this warranty policy.

If the packaging has not been opened, UNICSKIN will accept the refund within 14 days of the date of purchase at the customer's expense and risk. The device must be returned in the same packaging with the original proof of purchase. Please contact Customer Service for details of the returns address by emailing customerservices@unicskin.com, or by calling 00 34 912 797 377.

In the event of a manufacturing fault during the first two (2) years from the date of purchase, the return of any of our Techno-Beauty products may result in a like-for-like replacement (in which case the cost of return shall be paid by UNICSKIN), a refund or a reduction in the product price if the fault detected is minor or if the value of the goods at the time of making the claim does not warrant a replacement due to the time that has passed. In any case, it will be necessary to keep and present the product receipt or invoice with the date of purchase clearly visible.

Furthermore, UNICSKIN reserves the right to carry out the appropriate checks on the product purchased and may request photographs or evidence of the alleged defect, or that the faulty product is returned, with the shipping costs initially paid by the customer, and will reimburse them once UNICSKIN has verified that the fault detected is indeed a manufacturing defect.

Products damaged due to inappropriate storage conditions, poor maintenance, or any other type of misuse by the customer are excluded from this warranty policy. The product complies with the legal and regulatory standards applicable in the European Community at the date of manufacture and distribution.

UNICTECH EYE MASK ACCESSORIES

UNICTECH EYE MASK is lightweight, easy to use and convenient to take anywhere. Your device includes different accessories to make your life easier.



User
Manual



USB Charger



Case



Pouch



Cleaning wipe

COMPLEMENTARY PRODUCTS

Complete your eye contour routine with our dermocosmetic products designed to enhance the results.



UNICEYES TRIPLE ACTION

Latest generation expert treatment with 93% natural ingredients and a triple action to correct and prevent the signs of ageing: anti-wrinkle, anti-dark circles and anti-puffiness. Also decongests and soothes inflammation for a rested appearance. Brighter, more youthful eyes.

Cool effect serum texture for greater results.

100% NON-IRRITANT - 100% OPHTHALMOLOGICALLY TESTED



UNICCOLLAGEN FLASH EYE MASK

Collagen patches with anti-ageing properties that reduce wrinkles, dark circles and puffiness.

Rested, bright and beautiful eyes in an instant.

CUSTOMER SERVICE

If you have any questions or comments you would like to share with us, please email us at: customerservices@unicskin.com.

UNICSKIN

YOUR EFFECTIVE SOLUTION

GLOBAL EXPERTS IN TECHNO-BEAUTY & HOME BEAUTY
EFFECTIVE RESULTS PROVEN BY SCIENCE
SUSTAINABILITY & +94% NATURAL ORIGIN FORMULAS
INNOVATIVE LUXURY AT YOUR FINGERTIPS
